

## SMART DOG PARKING FACILITY

119-065



### PRODUCT DESCRIPTION

**A lead restraining hook and mounting plate, which may be securely attached to any suitable surface. A simple way to tether a dog**

Ideal for veterinary practice, kennels or groomers the Smart Dog Parking Facility can be easily fitted in receptions, waiting rooms, wards, prep rooms, exercise runs and consulting rooms. With no spring-loaded fastenings to contend with, the hand loop of a lead is passed over the hook and restrained by the 'cow horn' cross piece. Dogs don't have to be taken off the lead to attach so the park can be used single-handed. In the reception or waiting room, the Dog Park allows owners to tether dogs safely while dealing with payments or enquiries – it can be difficult to do this and control an excited dog at the same time! All metal parts are stainless steel.

The mounting plate is made from ABS plastic making the unit durable and damage resistant. Each unit is supplied complete with screws, wall plugs and fitting instructions. The Dog Parking Facility should not be used unattended in public areas. Apart from the risk of theft, there is the possibility that someone might frighten or untie the dog resulting in an accident.

#### Features

- Secure means to tether a dog
- Simple to fit and use
- Maintenance free

### ADDITIONAL INFORMATION

Animal

Small Animal

SKU	119-065
Overall Product Dimensions (mm)	110w x 130h
Colour/Style	White
Fit	N/A

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights