

## 1755MM LONG 90° BRACKET FOR KENNEL DOOR

131-119



### PRODUCT DESCRIPTION

## 90° bracket for kennel door to secure to partition wall above coving

- 1755mm Long

### ADDITIONAL INFORMATION

<b>Animal</b>	N/A
<b>SKU</b>	131-119
<b>Colour/Style</b>	Stainless Steel
<b>Brand</b>	Burtons
<b>Fit</b>	Walk-In (Glass), Walk-In (Mesh)
<b>Type</b>	Walk-in Accessories, Walk-in Door

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights