

## A4 CAGE CLIP BOARD

165-400-1290



### PRODUCT DESCRIPTION

**A4 Portrait Stainless steel clipboard that is designed to last a life-time, securely holding patient information to your cage doors**

- Made from a durable & hygienic T504 stainless steel allowing the clipboards to be easily wiped down and disinfected.
- Product includes bracket to fit a standard mesh cage and kennel door.
- Designed for Burtons cages

### ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	165-400-1290
<b>Colour/Style</b>	Stainless Steel
<b>Brand</b>	Burtons
<b>Fit</b>	Cage
<b>Type</b>	Bank Accessories

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights