

CHISON SONOBOOK 6 VET ULTRASOUND UNIT

USCH-SONOBOOK6-BUNDLE



Category: Tour Products > Ultrasound > Diagnostic Imaging

Price: £0.00

Date Printed: 29/03/24

PRODUCT DESCRIPTION

Designed to provide the best solution for various challenging veterinary environments. The SonoBook 6 features a compact, robust laptop design, weighing just 5.5kg and easily folding into a case for easy portability

With multiple advanced technologies, a 15 inch HD LED monitor with rotatable screen and 2hrs battery life in active mode makes this top of the range ultrasound machine incredibly good value for money, whilst still providing you with detailed, efficient examination using the latest technologies.

Key Features:

- 15 inch HD LED monitor
- Rotatable screen – 120 degrees
- Compact and robust laptop design – Metal casing, solid yet light
- Just 5.5kg for easy mobility
- Slim design – 70mm thick
- 2hrs battery life in active mode for outdoor Vet scanning – indication light – Easy to change
- Rich clinical application
- Professional Vet software
- Intuitive workflow
- Docking: DVI, Foot-switch, Video-out, Remote, S-Video

Advanced Technologies:

FHI

- Uses different transmissions and receiving methods for different sized animals – Maximising the resolution without losing the penetration.

X-Contrast

- Enhance the contrast resolution of different tissues by improving the signal-to-noise ratio
- Automatically detect different tissues and improve the contrast accordingly

Q-Image

- Innovative algorithms triple the image enhancement results to a quantum level
- Images are sharper and more vivid

Q-Beam

- Different to the traditional dual phased digital former, the quad beam is used to improve the processing speed and frame rate.

Q-Flow

- Maximise the colour sensitivity by automatically distinguishing different tissues and detecting colour flow accordingly.
- Better, colour sensitivity even for low-velocity flow and small vessels.

Professional vet software

- Body mark
- Comment library
- Calculation package

Available accessories to complete your system:

Dedicated trolley available with a small footprint – Includes anti-theft lock, accessory basket, printer socket, pneumatic height adjustment, sonotriple connector, ergonomic hand rest.

Accessories available: PC printers, video printers, dedicated carrying case, waterproof keyboard case

PRODUCT OPTIONS & ADD ONS

Ultrasound System

1x SonoBook 6 Vet Ultrasound Machine

Add your probes

1x L12-V Linear Probe
1x C3-V Convex Probe
1x P2-V Phased Array Probe
1x R7-V Transrectal Probe
1x MC3-V
1x L8M5-V Linear Probe
1x MC6-V Micro-Convex Probe

ADDITIONAL INFORMATION

Animal	Mixed
SKU	USCH-SONOBOOK6-bundle
Brand	Chison
Ultrasound screen size	15 inch
Ultrasound colour	Colour

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights