

**TERUMO NEEDLE BOX OF 100**

400-0-NEEDLES

**PRODUCT DESCRIPTION**

**Terumo's needles use a unique double bevel design and precision manufacturing to provide ultra sharp needles, which translate into greater comfort for the patient**

This product contains no components made of natural rubber latex and their tamper-evident blister packs result in the minimal waste disposal

- Terumo Needle, supplied in a box of 100
- UltraSharp Needles
- Tamper Evident Peel Packs

**ADDITIONAL INFORMATION**

<b>Animal</b>	N/A
<b>SKU</b>	400-0-needles

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights