

OPEN PALM MITTENS

146-854



PRODUCT DESCRIPTION

Open Palm Mitten with 0.50Pb Lead - Perfect for X-ray

Key Features:

- Lead equivalence of 0.50mm
- Core material of 'Lightweight Lead'
- Please refer to your Radiation Protection Advisor for advice on your particular usage
- Manufactured using multiple thin layers of protective sheet for maximum flexibility and suppleness
- Navy nylon covering materials
- New fully bound comfort binding
- Item comes as a pair

ADDITIONAL INFORMATION

| | |
|---------------|---------|
| SKU | 146-854 |
| Animal | N/A |

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

