

X-R STATIC UNIVET PACKAGES

146-620-30K-PACKAGES



Category: DR Systems > Radiography Systems > CR Systems

PRODUCT DESCRIPTION

The X-R Static UNIVET 37-400 is a cost-effective versatile high-frequency X-ray system that delivers high-quality radiographs

Thanks to the floating table top with electronic brake and integrated pre-selection of 352 standard organ programs, handling X-ray system is kept convenient and user-friendly. The X-R Static UNIVET 37-400 is available as standalone with a tablet for the species/organ selection and kV/mAs setting, also available as a system solution with an EXAMION cr system or dr system.

Key Features

- Fast and precise X-raying
- Easy installation and wall-plug suitable
- 352 standard organ programs
- Minimal Radiation Exposure
- Floating tabletop with electromagnetic brakes

Choose a CR or DR system to complete your purchase

Specifications

High-frequency X-ray device	
Output power	37.5 kW / 400 mA
Voltage range	40 - 100 kV
Current range	0.1 - 100 mAs
Switching time	0.001 - 1.25 s
Mains connection	230 V - 50 Hz 1-phase (wall plug)
Capacitor technology	3 kW input, 37.5 kW output
Programs	352 pre-programmed anatomical stations
Foot release switch	2-stage with safety function
Tabletop	Floating two movements
Table height	85 cm
Tabletop surface	160 cm x 68 cm, offset ± 28 cm
Tube	X-ray tube with rotary anode and dual focus 0.8 / 1.8 mm
Collimator	Continuously adjustable light beam type with cross indicator; manual rotation, halogen lamp, timer
Optional	X-R Static Univet as Standalone with a Tablet for the species/organ selection and kV/mAs Setting.
Available Options	<p>X-CR Static Univet with Tablet and an EXAMION CR System: X-CR Tabletop Vita, X-CR Tabletop Smart, Colenta HighCap Xr.</p> <p>X-CR Static Univet with an EXAMION CR System: X-CR Tabletop Vita, X-CR Tabletop Smart, Colenta HighCap Xr.</p> <p>X-DR Static Univet available with wired detectors: X-DR XL PE / X-DR XL RA - 43 cm x 43 cm (17" x 17"). Other detectors on demand.</p>
Warranty	12-month warranty Extended warranty on request

PRODUCT OPTIONS & ADD ONS

X-CR Univet System

1x Examion X CR Tabletop Vita Flex
 1x Examion X-CR Smart
 1x Examion Highcap Xr CR System

X-DR Univet System

1x Examion X DR XL PE Wired Detector Panel

X-R Static

1x Examion X-R Static UNIVET 37-400 Combined Table Generator

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	146-620-30K-packages
Brand	Examion
Catagory Sub Heading	Build Your System
System Type	Floor Standing
Upgrade Panel	Yes

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights