


## ACRYLIC DOOR FOR BURTONS CAGE

165-480-\*BUNDLE



### PRODUCT DESCRIPTION

**Taking the place of your standard door the acrylic door allows for much easier patient monitoring and helps reduce the spread of disease**

#### Features

- Constructed from 6mm clear acrylic
- Built in vent holes
- Stainless steel hinge and latch mechanism
- Simple to fit with a full range of sizes available

### ADDITIONAL INFORMATION

|                             |                        |
|-----------------------------|------------------------|
| <b>Animal</b>               | Small Animal           |
| <b>SKU</b>                  | 165-480-*bundle        |
| <b>Colour/Style</b>         | Clear                  |
| <b>Brand</b>                | Burtons                |
| <b>Catagory Sub Heading</b> | Various Size Avaliable |
| <b>Type</b>                 | Bank Accessories, Cage |

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights