


ACRYLIC DOOR FOR BURTONS CAGE

165-480-*BUNDLE



PRODUCT DESCRIPTION

Taking the place of your standard door the acrylic door allows for much easier patient monitoring and helps reduce the spread of disease

Features

- Constructed from 6mm clear acrylic
- Built in vent holes
- Stainless steel hinge and latch mechanism
- Simple to fit with a full range of sizes available

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	165-480-*bundle
Colour/Style	Clear
Brand	Burtons
Catagory Sub Heading	Various Size Avaliable
Type	Bank Accessories, Cage

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights