

BOWL HOLDER RING WITH INTEGRAL MOUNTING BRACKET

165-420-23-BUNDLE



PRODUCT DESCRIPTION

Strong bowl holder ring with cage mounting bracket

(Bowl not included)

- Save space on the cage floor
- Bowl is easily added / removed
- Attached securely to mesh doors
- Open door and gain full access to bowl, away from patients
- Available to fit a range of Burtons bowls
- Strong steel finish

ADDITIONAL INFORMATION

Animal	Mixed
SKU	165-420-23-bundle
Colour/Style	Silver
Brand	Burtons
Catagory Sub Heading	Available In Various Sizes
Fit	Cage
Type	Bank Accessories

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights