

## SMART IV CATHETER GUARD

119-07-CONFIG



### PRODUCT DESCRIPTION

---

## Lightweight flexible plastic shield designed to cover the site of the intravenous catheter

#### Key Features

- Helps keep IV catheter sites safe
- Reduces need for bandaging
- Simple to use

Reducing the need for excessive taping and bandaging, our IV Catheter Guard is well tolerated and can help reduce the frequency of mega-paw swelling. Easy to apply, simply select the correct size guard allowing good coverage of the catheter insertion site. Loosen velcro strap and place guard centrally over the catheter site. Wrap velcro straps around the limb and fasten. Available in three sizes.

#### Available Sizes

Select correct size guard to allow good coverage

- Medium - 9cmL x 9cmW
- Large - 13cmL x 10cmW
- Extra Large - 17cmL x 12cmW

### ADDITIONAL INFORMATION

**Animal**

Small Animal

SKU

119-07-config

Colour/Style

Black

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights