

DERMAFREEZE FOAM TIPPED APPLICATOR REFILL PACK

128-31-REFIL



Category: PPE & Accessories > Consumables

Price: £20.53 - £41.06
Date Printed: 29/03/24

PRODUCT DESCRIPTION

Dermafreeze Cryosurgery foam tipped applicator refill packs

- Pack sizes x20
- 2mm and 5mm sizes

PRODUCT OPTIONS & ADD ONS

Select Size

- 1x 2mm Dermafreeze Foam Tipped Applicator Refill Pack
- 1x 5mm Dermafreeze Foam Tipped Applicator Refill Pack

ADDITIONAL INFORMATION

Animal	N/A
SKU	128-31-refil

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights