

BURTONS O₂ ANAESTHESIA MACHINE

102-012-BUNDLE



PRODUCT DESCRIPTION

Great value single gas mobile anaesthetic machine

Including a 0-15 l/min O₂ flowmeter and O₂ flush safety valve as standard, and complete with a 5 star mobile base. This anaesthetic machine offers exceptional value for money whilst still being built to last. Built to carry 1 x 'F' size cylinder it still remains compact and easy to move around the practice, whilst being easily stored when not in use.

Key Features:

- Compact Design
- Stable 5 Star Base
- Holds one F size bullnose cylinder
- 0-15 l/min O₂ Flowmeter
- Bullnose O₂ Regulator
- Emergency O₂ Flush
- Corner unit outlet (BS22mm Male)
- Overall size: 61 x 61 x 114cm WDH

Don't forget your Vaporiser

The Series 3 is a temperature compensating, out of circuit vaporiser in a traditional design, familiar to most veterinary practitioners. Calibrated specifically for the specified anaesthetic drug using a laser refractometer.

Fitting Type

Choose from the standard taper mounted fitting (cage mount) (UNI) or the quick release (QR) option; allowing you to easily swap the vaporiser to different anaesthetic machines throughout the practice. Additional quick release compatibility blocks will be required for

each anaesthesia machine. (Available separately)

Key filled or Funnel (Pour fill)

Key filled (KF) systems require the addition of key filler adaptors but are generally safer and have a lot less wastage compared to standard funnel filled (FF) vaporisers, saving you money in the long run.

Manufactured to last

We don't just supply it! It was designed and built by Burtons in the UK at our principal manufacturing facility in Maidstone, Kent Using the latest CNC's and laser manufacturing processes and high quality stainless steel materials this product has been engineered from the ground up to last!

Vaporiser Sold Separately

ADDITIONAL INFORMATION

| | |
|--|-----------------------------|
| Animal | Small Animal |
| SKU | 102-012-bundle |
| Overall Product Dimensions (mm) | 610mm W x 610mm D x1140mm H |
| Colour/Style | White |
| Gas Type | O ² |
| Gas Supply | 'F' Size Cylinder(s) |
| Flowmeter Type | Flowmeter |
| Concentrator | No |
| Optional O2 Alarm | No |
| Optional N2O | No |
| Emergency O2 Flush | Yes |

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

