


## ANAESTHESIA PATIENT TROLLEY

200-113-BUNDLE



## PRODUCT DESCRIPTION

### Increase mobility with the purpose built anaesthesia patient trolley complete with aga holder

Constructed from high quality, medical grade, stainless steel this combination trolley offers a strong, sturdy base for procedures with the added benefit of a built in anaesthetic machine, charcoal canister and 'F' size cylinder holder. Meaning you can keep your patient anaesthetised whilst on the move without the need for an additional trolley.

Additional benefits include door runners for high speed transfer of patients, a high accuracy rotameter and emergency O<sub>2</sub> flush as standard.

#### Key Features:

- Stainless Steel Frame
- Holds one "F" size cylinder
- Supplied with one regulator
- Rapid cylinder change
- High accuracy rotameter 10 l/min O<sub>2</sub>
- Emergency O<sub>2</sub> flush
- Safety pressure relief valve and non-return valve
- Built in aga holder
- Can also be used with piped gas systems
- Dimensions: 740mmW x 1560mmL x 1070mm H (see diagram of dimensions above)

### Don't forget your Vaporiser

The Series 3 is a temperature compensating, out of circuit vaporiser in a traditional design, familiar to most veterinary practitioners. Calibrated specifically for the specified anaesthetic drug using a laser refractometer.

### Fitting Type

Choose from the standard taper mounted fitting (cage mount) (UNI) or the quick release (QR) option; allowing you to easily swap the vaporiser to different anaesthetic machines throughout the practice. Additional quick release compatibility blocks will be required for each anaesthesia machine. (Available separately)

### Key filled or Funnel (Pour fill)

Key filled (KF) systems require the addition of key filler adaptors but are generally safer and have a lot less wastage compared to standard funnel filled (FF) vaporisers, saving you money in the long run.

### Manufactured to last

We don't just supply it! It was designed and built by Burtons in the UK at our principal manufacturing facility in Maidstone, Kent Using the latest CNC's and laser manufacturing processes and high quality stainless steel materials this product has been engineered from the ground up to last!

*Vaporiser Sold Separately*

## ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	200-113-bundle
<b>Overall Product Dimensions (mm)</b>	740mmW x 1560mmL x 1070mm H
<b>Colour/Style</b>	Stainless Steel
<b>Gas Type</b>	O <sup>2</sup>
<b>Gas Supply</b>	'F' Size Cylinder(s)
<b>Flowmeter Type</b>	Graduated Rotameter
<b>Concentrator</b>	No
<b>Optional O2 Alarm</b>	No
<b>Optional N2O</b>	No
<b>Emergency O2 Flush</b>	Yes

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights