

## WELCH ALLYN PNEUMATIC DESK DIAGNOSTIC SET

WA71784



## PRODUCT DESCRIPTION

Welch Allyn pneumatic NiCad rechargeable desk diagnostic set includes:

### **Coaxial Ophthalmoscope:**

- 6 apertures for all general and specialist use
- 28 lenses -25 to +40 diopters for better resolution
- Rubber brow rest prevents scratching of eyeglasses
- Illuminated lens dial clearly identifies settings
- Coaxial optical system for shadow-free spot and easier entry into undilated pupil
- Crossed linear polarizing filter to reduce corneal reflection
- Red-free filter for use with any aperture
- 18 possible aperture combinations for greater versatility
- Sealed optics for better dust protection

### **Pneumatic Otoscope:**

- Rotating lens with ultra-seal for pneumatic otoscopy
- Fiber optic illumination for cool light

### **3.5 V Rechargeable Handle:**

- Economical—long-lasting rechargeable nickel-cadmium battery included
- Sure grip—smooth, nonabrasive knurled finish ensures a no-slip grip
- Chrome-plated brass handles with durable, long-life construction
- Safety switch to prevent accidental draining of battery
- Metal head connectors for durability

### **3.5 V Universal Desk Charger:**

- Advanced circuitry for constant automatic safe charging and optimum lamp performance and life
- Compact for convenience and mobility
- Compact charger base measures only 3-3/8" (86 mm) x 7-1/4" (184 mm)
- Optional wall brackets for permanently mounting to wall or cabinet

The pneumatic NiCad rechargeable desk diagnostic set diagnostic Set also includes 3 x polypropylene specula and (4,5 & 7mm).

## ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	WA71784
<b>Generic Size</b>	Small
<b>Colour/Style</b>	Black
<b>Brand</b>	Welch Allyn
<b>Digital</b>	No
<b>Diagnostic Head Type</b>	Pneumatic Otoscope

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights