

## DE-IONISING RESIN - 1LTR FOR USE WITH REFILLABLE FILTER

139-000-1



### PRODUCT DESCRIPTION

---

DI Mixed Bed Resin is used in our refillable cartridges. Demineralisation Resins are used to reduce TDS (polishing the water).

The TDS from a RO unit is not always as low as required so using DI resin will bring the TDS down to almost zero

### ADDITIONAL INFORMATION

SKU	139-000-1
-----	-----------

#### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights