

## DE-IONISATION RESIN FILTER 10" (EMPTY)

139-000-1-1



### PRODUCT DESCRIPTION

This refillable deionisation filter is ideal for fitting on to any commercial RO system as a final polishing filter.

Refill this filter with our Semiconductor Grade Mixed bed ion Exchange Resin and ensure 0 TDS at your system outlet.

Service life is dependant on the quality of the water supplied to this filter. It is advisable to use an [inline TDS Meter](#) with this filter.

### ADDITIONAL INFORMATION

|                     |             |
|---------------------|-------------|
| <b>SKU</b>          | 139-000-1-1 |
| <b>Animal</b>       | N/A         |
| <b>Colour/Style</b> | White       |

#### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

