

**PORTABLE ELECTRIC MICROMOTOR KIT**

118-880

**PRODUCT DESCRIPTION****Portable motor with E type electric micromotor with Jack plug fitting**

- Motor contains re-chargeable Li-Ion battery which can operate for 4-6 hours on a single charge.
- Motor speed 0-30,000 RPM
- Power output 12 v / 25 Watt/
- Small and lightweight ideal for domiciliary use
- 12 months warranty.

**ADDITIONAL INFORMATION**

<b>Animal</b>	N/A
<b>SKU</b>	118-880
<b>Colour/Style</b>	White
<b>Examination Type</b>	Oral Examination, Small Furies, Tooth Extraction

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights