

HUGER AGVE-68 SERIES BRONCHOSCOPES

END-HGVG-BRONCHOSCOPES



PRODUCT DESCRIPTION

HUGER AGVE-68 series GI endoscopes employ advanced CMOS image sensor and combined with the compact VIS-68 video Huger endoscopy systems (with superior algorithm technology and refined functions) present exceptionally high-resolution images with clarity in every detail.

- 1240*720p CMOS high resolution endoscope
- New generation handle, smaller, lighter and of remote control button and enhanced maneuverability
- Smoother, flatter and waterproof endoscope body for easier cleaning and disinfection

Model	68PL	68Q			
Field of view	120°	120°			
Direction of view	Forward viewing	Forward viewing			
Depth of view	3~100mm	3~100mm			
Distal end diameter	3.8mm	5.6mm			
Insertion tube diameter	3.8mm	5.6mm			
Instrument channel diameter	No channel	2.0mm			

Angulation range	U/D210°/90°;L/R100°	U/D210°/90°;L/R100°			
Working length	600mm	600mm			

ADDITIONAL INFORMATION

Animal	Mixed
SKU	END-HGVG-BRONCHOSCOPES
Brand	Huger
Procedure Type	Bronchoscopy, Cystoscopy, Respiratory Examination, Rhinoscopy
Biopsy Channel	1.4mm

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights