

## MEDICAPTURE MVR VIDEO ARCHIVING UNIT (HD)

END-MCMVR409



### PRODUCT DESCRIPTION

---

**Based on the exclusive iMave pro platform, the MVR Pro is the most advanced, premium-quality medical video recorder with full patient workflow capability and an easy built-in touch screen designed to benefit today's modern operating rooms and surgeries.**

#### Key Features

- Next generation video recorder – Allowing for multi input recording in crystal clear Full HD 1080p 60 resolution.
- Data protection – due to its data life management
- Full patient workflow – Including DICOM option for easy and fast retrieval of patient information, transfer of video files and network storage
- Simple to use – The unique built-in touchscreen includes intuitive icons making the recorder very easy to use.
- Full network connectivity
- Multi-input recording
- Multi-Storage/ DLM
- Optimized recording
- DICOM App
- Streaming App
- IT Remote App
- Edit and print App

### ADDITIONAL INFORMATION

<b>SKU</b>	END-MCMVR409
<b>Animal</b>	Mixed

**Generic Size**

Small

**Colour/Style**

White

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights