

## ISOFLURANE KEY FILLER BOTTLE ADAPTOR

102-921



### PRODUCT DESCRIPTION

## High Quality Key fill adapter facilitating the filling of key fill style vaporisers safely whilst minimising waste

Each adapter is agent specific and colour coded for the appropriate anaesthetic agent.

- Purple = Isoflurane.
- Red = Halothane
- Yellow = Sevoflurane
- Requires Key Fill Style Vaporiser

### ADDITIONAL INFORMATION

**Animal**  
**SKU**  
**anaesthesia\_vap\_anaesthetic**

Small Animal

102-921

Isoflurane

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights