

ISOFLURANE KEY FILLER BOTTLE ADAPTOR

102-921



Category: Anaesthesia Consumables > Anaesthesia Consumables > Vaporisers > Anaesthesia

Price: £56.23
Date Printed: 17/04/24

PRODUCT DESCRIPTION

High Quality Key fill adapter facilitating the filling of key fill style vaporisers safely whilst minimising waste

Each adapter is agent specific and colour coded for the appropriate anaesthetic agent.

- Purple = Isoflurane.
- Red = Halothane
- Yellow = Sevoflurane
- Requires Key Fill Style Vaporiser

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	102-921
Anaesthetic	Isoflurane

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights