

## MICRO SD CARD FOR LABCOLD REFRIGERATORS

125-047



### PRODUCT DESCRIPTION

To get the most out of your touchscreen, Labcold Fridge, you will need to purchase a Labcold specific Micro SD card. This will allow you to record all data and transfer it to your computer

- Minute by minute, 24/7, 365 days a year temperature data for your records
- Collects data from both temperature probes

The IntelliCold® controller in your pharmacy refrigerator will automatically store readings from the two temperature probes plus any alarms. It will also store information about the fridge itself, such as whether the door was open and even if the fan or compressor are running and, if so, for how long. This data is collected and saved every minute, 24/7, 365 days of the year. To view this information as a csv on your computer, you must have an encoded Micro SD card from Labcold to transfer it from the fridge to your desktop.

As long as the refrigerator has power, this data will be recorded and stored, giving you all the information you need for regulatory inspections and meticulous and comprehensive records. In addition, should you experience an issue you can check, minute by minute, whether the temperature inside the refrigerator was affected.

Designed by Labcold for our IntelliCold® Pharmacy Range RLDF/GXX19, the micro SD card comes complete with an adapter, so is suitable for a computer which only has an SD card slot. There is also a case to keep the Micro SD and adapter together. To ensure utmost security, only micro SD cards from Labcold will work with your fridge and the first time the card is used, it will lock to that fridge to prevent data from being lost or recorded from other refrigerators.

**\*If you already have your Labcold fridge, please supply us with your fridge serial number when ordering\***

### ADDITIONAL INFORMATION

<b>SKU</b>	125-047
<b>Brand</b>	Labcold

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights