

MODIFIED AYRES 'T' COMPLETE CIRCUIT

103-050

**PRODUCT DESCRIPTION****A modified ayres T non-rebreathing circuit typically for use with patients of below <10kg body weight**

- Suggested minimum flow rate 2.5 times minute volume
- 1.2mtr feed tube
- 22M/15F patient connector
- Supplied with 30mm scavenge bag adaptor, 0.5 litre bag
- Typically for patients < 10Kg
- Safety factor 2.5 x minute volume
- Mapleson F classification

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	103-050
Durability	Non-Disposable
Tubing Length	1.2 Meter
MRI Compatible	No
Circuit Type	Modified Ayres T

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights