

NON DISPOSABLE CO-AXIAL BAIN CIRCUIT

103-055



PRODUCT DESCRIPTION

A non rebreathing circuit with co-axial tubing; fresh gas is supplied through the centre bore whilst exhaust gases are exhaled through the outer tubing

Recommended for use with patients above 10Kg body weight.

- Suggested minimum flow rate 2.5 times minute volume
- 1.6 metre Co-axial tubing set
- 2M/15F patient connector
- 30mm scavenge valve
- 2 litre rebreathing bag
- Typically for patients > 10Kg
- Safety factor 2.5 x minute volume
- Mapleson D Classification

ADDITIONAL INFORMATION

| | |
|-----------------------|----------------|
| Animal | Small Animal |
| SKU | 103-055 |
| Durability | Non-Disposable |
| Tubing Length | 1.6 Meter |
| MRI Compatible | No |
| Circuit Type | Bain |

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights