

CYCLO-FLO CIRCLE CIRCUIT

103-080



Price: £656.00

Date Printed: 16/04/24

Category: Product Imports > Anaesthetic Circuits > Anaesthesia

PRODUCT DESCRIPTION

Designed specifically to provide the lowest resistance possible

Circles have become the standard circuit of choice for Vet's so we have designed the Cyclo-Flo to provide reliability, economy and low resistance.

The Cyclo-Flo Circle Breathing System is supplied with a 2-litre rebreathing bag, 22mm x 1.6m parallel patient tubing, V mounting clamp bracket and operator manual.

Key Features

- Designed specifically to provide the lowest resistance possible
- Compatible with virtually any anaesthetic machine
- Large capacity (1.3kg) easy fit canister
- Canister septum reduces absorbent wastage
- High visibility one-way valves
- Simple and minimal maintenance
- Medical V type mounting - easy to fit and remove
- Suitable for bodyweight range 7-150kg
- MRI compatible version available

ADDITIONAL INFORMATION

Animal	Mixed
SKU	103-080
Durability	Non-Disposable

MRI Compatible

No

Circuit Type

Circle System

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights