

## DISPOSABLE MODIFIED AYRES T CIRCUIT WITH GAS SCAVENGE ADAPTOR

103-210



### PRODUCT DESCRIPTION

## Disposable Modified Ayres T Circuit With Gas Scavenging Adaptor For Bag

- High quality
- 1.8mtr feed tube
- 22M/15F patient connector
- Supplied with 30mm scavenge bag adaptor
- 0.5 litre bag
- Typically for patients < 10Kg.
- Safety factor 2.5 x minute volume

### ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	103-210
<b>Durability</b>	Disposable
<b>Tubing Length</b>	1.8 meter
<b>MRI Compatible</b>	No
<b>Circuit Type</b>	Ayres T

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights