

**TOFFELN FLEXIKLOGS**

0723-FLEXIKLOGS-CONFIG\*

**PRODUCT DESCRIPTION****Supports good posture**

FlexiKlog is designed for comfort. It has a built-in arch and a slightly raised heel to help keep your body in the correct position, which helps you stay more comfortable all day long. There's also plenty of room for your toes to spread out, thanks to the high toe box.

- Soft Leather Insole - cushioned for ultimate comfort.
- Anti-Static Sole - Shock absorbing sole specifically incorporated to meet static sensitive situations.
- Ergonomically Designed - fully flexible sole is designed to provide day-long comfort.
- Leather Upper - PU coated easy-to-clean uppers.
- High Toe Box - gives room for your toes to spread

**ADDITIONAL INFORMATION**

<b>Animal</b>	N/A
<b>SKU</b>	0723-flexiklogs-config*
<b>Male/Female</b>	Unisex
<b>Brand</b>	Toffeln

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights