


LITTMANN LIGHTWEIGHT NURSE STETHOSCOPE

117-671-CONFIG1*

**PRODUCT DESCRIPTION**

The 3M™ Littmann® Classic II S.E. Stethoscope is a dependable, hard working clinical tool for physical assessment and monitoring of your patients. It's the choice of millions of care providers who are dedicated to achieving the best patient outcomes.

Key Features

- Versatile chestpiece design provides high acoustic sensitivity with both a tunable diaphragm and open bell.
- Tunable diaphragm: Hear high or low frequency sounds by slightly adjusting pressure on the chestpiece.
- Non-chill rim and diaphragm provide patient comfort.
- Soft-sealing eartips provide an excellent acoustic seal and comfortable fit.
- Headset is easily adjusted for individual fit and comfort. Angled eartubes align with ear canals.

ADDITIONAL INFORMATION**Animal**

Small Animal

SKU

117-671-config1*

Stethoscope Type

Lightweight SE

Chestpiece Style

Standard

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights