

## SOLITAIRE ANAESTHESIA MACHINE - O<sup>2</sup> ONLY CONCENTRATOR READY - CLEARANCE

127-100-CON-CL



### PRODUCT DESCRIPTION

---

## Wall mounted anaesthetic machine from the Burtons Value Range

A Low cost alternative to other anaesthetic machines the solitaire is made from durable powder coated mild steel with high quality built in 0-15lpm rotameter and emergency O<sup>2</sup> flush. The slim-line mobile base with stable 5 star design enables the machine to be easily transported around the practice, whilst being compact enough to be out the way during storage.

Add an optional E or F size cylinder kit for easy cylinder storage on-board the machine.

#### Features

- Monitor Shelf
- Stable 5 star mobile base
- 0-15lpm O<sup>2</sup> Rotameter with emergency O<sup>2</sup> Flush
- Corner unit outlet ( BS22mm Male)
- Compact Design
- Available in piped gas or E or F Size Cylinder Versions
- Overall size: 122 H x 62cm Diameter
- Designed, developed and manufactured by Burtons in the UK

- Pipeline Lead & Vaporiser Sold Separately
- Machine & concentrator ONLY

Clearance products may be non-pristine, for example:

- On occasion some items may show signs of minimal use, which may include minor cosmetic marks which do not affect the functionality in any way.
- Although most of our clearance products are sold in their original manufacturers packaging, there may be occasions, due to packaging damage, where plain white or brown packaging may be used

Clearance units sold with 6 month standard guarantee, unless otherwise specified. This does not affect your statutory rights.

Please see our terms and conditions for more information.

## ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	127-100-CON-cl
<b>Colour/Style</b>	White
<b>Gas Type</b>	O <sup>2</sup>
<b>Gas Supply</b>	'E' Size Cylinder (s), 'F' Size Cylinder(s), Piped Gas
<b>Flowmeter Type</b>	Graduated Rotameter
<b>Concentrator</b>	No
<b>Optional O2 Alarm</b>	No
<b>Optional N2O</b>	No
<b>Emergency O2 Flush</b>	Yes

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

