

BURTONS PRO-LIFT MONO TABLE BASE - CLEARANCE

225-207-CL



PRODUCT DESCRIPTION

The latest addition to our range of tables, this ultra-quiet, electric examination table excels in functionality, stability and longevity

Get your patient to the correct position with ease with the Burtons Pro-Lift mono. With adjustable height controlled using the electric powered ram, reducing repetitive strain and ensuring maximum comfort. Its simple controls enable you to adjust the table height from 690mm to 1090mm and the built-in manual tilt allows you to get optimum angle and positioning for all examinations.

With a standardised table top fitting this table can be configured for your needs allowing for most of the Burtons table top range to be fitted with additional heated and lead lined options available. Simply select the top you need.

Electric Height Control

Enjoy ease of use around your practice with built-in battery, allowing you to work wherever you need, and don't be afraid of it losing power, the table can be used directly from mains as a backup.

The Frame

Made up of a single, heavy duty, electric ram the simple design of this table base allows for quick and easy cleaning. All supported on a heavy duty stainless steel base with four anti-static medical grade, lockable castors, ensuring smooth and easy manoeuvrability around the practice.

Manufactured to last

We don't just supply it! It was designed and built by Burtons in the UK at our principle manufacturing facility in Maidstone, Kent Using the latest CNC's and laser manufacturing processes and high quality stainless steel materials this product has been engineered from the ground up to last!

Key features -

- Electric Lift
- Hand Controls
- Battery As standard
- Manual Tilt
- Hygienic and easy to clean
- Ultra-quiet operation
- Easy to use controls
- Heavy-duty, high-quality stainless steel base.
- Twin wheels - lockable, anti-static castors.
- Minimum height: 47cm
- Maximum height: 95cm
- Maximum load: 200Kg

Clearance products may be non-pristine, for example:

- On occasion some items may show signs of minimal use, which may include minor cosmetic marks which do not affect the functionality in any way.
- Although most of our clearance products are sold in their original manufacturers packaging, there may be occasions, due to packaging damage, where plain white or brown packaging may be used

Clearance units sold with 6 month standard guarantee, unless otherwise specified. This does not affect your statutory rights.

Please see our terms and conditions for more information.

S/N: 055274

ADDITIONAL INFORMATION

| | |
|---------------------|-----------------|
| Animal | Small Animal |
| SKU | 225-207-cl |
| Colour/Style | Stainless Steel |
| Brand | Burtons |
| Mobile | Yes |
| Control Type | Handheld |
| Lift Type | Electric |

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights