

## ADAPTABLE2 ANAESTHESIA MACHINE

102-060-BUNDLE\*



## PRODUCT DESCRIPTION

### The same industry-leading Adaptable... in a new, smarter material!

**The Adaptable2 is manufactured in a solid uPVC material frame that is impervious to liquid and resistant to bacteria, making for a completely sealed, easy-wipe and waterproof anaesthesia machine that stands up to rigorous cleaning regimes and reduces contamination risk.**

Built upon the original Burtons AAA (Adaptable Anaesthesia Apparatus) frame concept, the Adaptable2 Anaesthesia System is designed for versatility and can form the basis of a variety of anaesthesia circuits. Perfect for piped gas, this machine can also accommodate two F-size oxygen cylinders in hinged clamps on the rear and an on-board oxygen concentrator on the base shelf. Any and all of these oxygen supply methods can be used without the need to reconfigure your machine's connectors.

For peace of mind in theatre, the Adaptable2 includes a number of carefully considered safety features. In cases of power outage, the manual changeover valve allows you to quickly and seamlessly switch from your on-board concentrator to your backup cylinders. The emergency flush button contains two integral safety valves for pressure relief and to prevent expired gases flowing back through the circuit to the vaporiser and beyond.

#### Key Features

- Frame made from a waterproof, antibacterial, solid uPVC - a continuous easy-wipe surface;
- High-accuracy 0-10l/min logarithmic rotameter;
- Emergency O<sub>2</sub> flush includes integral pressure relief and non-return safety valves;

- Connects seamlessly with the piped oxygen in your practice, F-size cylinders, or with an on-board concentrator;
- Rapid manual change-over valve for switching between concentrator and cylinder oxygen supplies – vital for power outages;
- Includes a monitor shelf, situated on the top of the frame;
- Can accommodate up to two vaporisers;
- Fully mobile, solid four-wheel base with independent wheel locks;
- Full warranty, supported by Burtons’ nationwide Service Team;
- Designed and manufactured by Burtons in the UK
- This machine frame is 1380mm High x 480mm Wide x 650mm Deep. With the optional Cyclo-Flo & Air Brake fitted the width increases to 630mm (See dimensions diagram above)
- F-size cylinder rings and regulator kit, concentrator and vaporiser(s) sold separately.

## ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	102-060-bundle*
<b>Gas Type</b>	O <sup>2</sup>
<b>Gas Supply</b>	'F' Size Cylinder(s)
<b>Concentrator</b>	No
<b>Optional O2 Alarm</b>	No
<b>Optional N2O</b>	No
<b>Emergency O2 Flush</b>	No

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights