

BURTONS ANAESTHETIC GAS ABSORBER - BOX OF 6

133-090



Category: Anaesthesia Consumables > Consumables > Anaesthesia Consumables > Anaesthesia Accessories > Anaesthesia > Scavenging

Price: £105.00
Date Printed: 06/06/23

PRODUCT DESCRIPTION**Box of 6 activated charcoal canisters for the absorption of anaesthetic gases**

The Burton's AGA is designed for the absorption of organic anaesthetic gases e.g. halothane and isoflurane. This method is not suitable for use with nitrous oxide because it is not absorbed by the charcoal. The weight of the canister must be monitored so that it can be replaced before it becomes saturated, once saturated any waste gases will simply be exhausted into the working environment.

Key Features

For absorption of gases such as halothane & isoflurane

Connected directly to the scavenging tubing from the patient circuit expiratory (scavenge) valve

Inlet connection: 22MM male

ADDITIONAL INFORMATION

Animal	N/A
SKU	133-090

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights