

## LARGE ANIMAL FOAL CIRCUIT

V7138



### PRODUCT DESCRIPTION

## Large Animal Foal Circuit

The Foal Circuit Set is commonly used with the LDS 3000 Large Animal Anaesthesia Machine. The hoses included in the set can also be used with SurgiVet brand small animal anaesthesia machines.

#### Kit Includes:

- One 3-litre bag
- Two 60-inch 22 mm hoses
- Three large animal endotracheal tube funnels
- Three funnel valve adapters
- One "Y" piece with Luer Fitting

### ADDITIONAL INFORMATION

<b>Animal</b>	Equine
<b>SKU</b>	V7138
<b>Durability</b>	Disposable
<b>Tubing Length</b>	60"
<b>MRI Compatible</b>	No

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights