

**OXYGEN DEMAND VALVE L063**

103-460-CONFIG\*

**PRODUCT DESCRIPTION****Designed to provide 100% oxygen on demand or by button activated supply**

This Resuscitation Demand Valve should be used in conjunction with a regulated oxygen supply.

- 22mm patient connection port
- 1/4 male gas inlet port
- 160 Litres per minute flowrate

**ADDITIONAL INFORMATION**

<b>Animal</b>	Small Animal
<b>SKU</b>	103-460-config*
<b>MRI Compatible</b>	No

**Good to know...**

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

