

## ANNEQUIN FOR ANAESTHESIA TRAINING

VSL-VS0800



### PRODUCT DESCRIPTION

**Annequin (animal mannequin) is a new dummy animal capable of accurately emulating the behaviour and responses seen in ventilated animals with respect to CO<sub>2</sub> elimination.**

The Annequin kit can be used with any standard ventilator and any form of capnograph. Annequin mimics a normal animal by having a constant rate of CO<sub>2</sub> production that must be removed by ventilation. However, you will need to supply the Annequin with a feed of 100% CO<sub>2</sub>. This is usually readily available in practice as this gas is used for insufflation during laparoscopies. Our kit includes a precision low flow (0-100ml) flow meter to allow you to accurately set the CO<sub>2</sub> flow into your the annequin.

With Annequin you can easily simulate hyper and hypocapnic states, investigate the effects of over and under-ventilation, as well as the effects of dead space and flow rates in breathing circuits.

This allows safe practical training for vets and nurses in responses to emergencies in ventilation.

Supplied annequins have a typical tidal volume of 150-200 mls and a compliance of approximately 15mls/cm.

Watch our video giving an Introduction to the Annequin Training Model [here](#)

#### The Annequin kit (VS0800) consists of the following:

- 1 x single Annequin fitted with a 'lung' , oral 15mm ET connector and CO<sub>2</sub> feed portal
- 2732 - 1 x 2.45m CO<sub>2</sub> feed line with male-male luer connections
- VS0750 - 1 x 0-100ml/min CO<sub>2</sub> flow meter on a stand
- VS0760 - 1 x luer-luer stop-cock
- VS0765 - 1 x 2m length 4mm polyethylene tubing with 1 x male luer connector

### ADDITIONAL INFORMATION

SKU

VSL-VS0800

<b>Brand</b>	Vetronic Services Ltd
<b>Animal</b>	N/A

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights