

BURTONS PIPELINE LEAD

102-778-BUNDLE

**PRODUCT DESCRIPTION****Burtons O² and N²O Pipeline Lead**

Suitable for use with Burtons anaesthetic machines or with pin index pipeline blocks for most other anaesthetic machines

- Various lengths available
- High quality connectors

ADDITIONAL INFORMATION

| | |
|-----------------------------|--------------------------|
| Animal | Mixed |
| SKU | 102-778-bundle |
| Catagory Sub Heading | Range of Sizes Available |
| MRI Compatible | No |

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights