

## GAS ANALYSER WASTE KIT

VSL-VS0058



### PRODUCT DESCRIPTION

The VS0058 Gas Analyser Waste Kit is comprised of:

- **2732** gas monitoring line, 1.2mm ID, male/male luer lock **x1**
- **1963** Straight connector 22M - 6mm OD stem and cap - 22F **x1**
- **1965** Straight connector 22M - 6mm OD stem and cap - 22M **x1**

### ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	VSL-VS0058
<b>Type</b>	Other
<b>Monitor Fit</b>	Vetronic Impact, Vetronic Lightning

#### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

