


F-SIZE CYLINDER KIT FOR MOBILE SOLITAIRE

127-150



PRODUCT DESCRIPTION

F Size Cylinder Kit

Fitting kit for adding 1 'F' size cylinder to your mobile solitaire anaesthetic machine.

Kit includes

- 2 x clamps
- 1 x regulator
- 1 x pipeline lead

ADDITIONAL INFORMATION

Animal	Mixed
SKU	127-150
Gas Supply	'F' Size Cylinder(s)
Concentrator	No
Optional O2 Alarm	No
Optional N2O	No
Emergency O2 Flush	No

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights