

IMEC 8 WITH PRINTER - ECG/TEMP/SP02/NIBP/RESP/CO2 - CLEARANCE

IMEC8VET-CO2-P-CL



PRODUCT DESCRIPTION

The iMec 8 Veterinary Monitor is a small, lightweight monitor that has many features packed into it, the perfect addition to any practice.

This neatly packaged monitor has a 8.4" High resolution touchscreen, weighs less than 3.2Kg (making it easily portable) and has up to 4 hours continuous monitoring and 120 hours trend / 48 hours waveforms reviewing.

Key Features

- 8.4" High-resolution Touch screen LED display with 8 waveforms
- Portable, less than 3.2 kg
- No-fan design for a quiet care environment
- Easy to use for intuitive operation
- Up to 4 hours continuous monitoring with large capacity Li-ion battery
- Parameters including CO2/ECG/Temp/SP02/NIBP/Respiration
- 360-degree visible alarm indicator
- 120 hours trend and 48 hours waveforms reviewing



ADDITIONAL INFORMATION

Animal

Mixed

SKU	IMEC8VET-CO2-P-cl
Colour/Style	White
Brand	Mindray
Monitor Parameters	CO ₂ , ECG, Multi-Parameter, NIBP, Respiration , SpO ₂ , Temperature
Screen Type	Touchscreen
Power Source	Mains, Rechargeable
Optional Modules	No
SD Card	No

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights