

**FROSTING TO BOTTOM HALF OF DOOR - 79.0CM W X 196CM H  
KENNEL DOOR - CLEARANCE**

131-101-G-F-L-CL

**PRODUCT DESCRIPTION**

Glass Covered Frosting to Bottom Half of Door on Outside - 79.0cm W x 196cm H Kennel Door H Fixing Brackets Not Included - LH Hinge (not reversible)

**S/N: 0000858967****ADDITIONAL INFORMATION**

<b>Animal</b>	Small Animal
<b>SKU</b>	131-101-g-f-l-cl
<b>Colour/Style</b>	Stainless Steel
<b>Brand</b>	Burtons
<b>Complete Bank</b>	No
<b>Material</b>	Stainless Steel
<b>Mobile</b>	No
<b>Door Type</b>	Half Glass / Half Mesh

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights