

**QS5 - ULTRAWAVE CLEANING BATH - 4.5L**

116-617

**PRODUCT DESCRIPTION**

## Ultrawave's Q-series baths allow accurate cycle parameter control and effective cleaning performance for all applications

The QS5 is a 4.5Litre heated, digital ultrasonic bath. The menu driven control panel allows simple adjustment of cycle time (0 to 99 minutes), temperature (ambient to 70°C), ultrasonic power and frequency mode for precise and accurate processing.

Combining features such as degas mode, fluid level sensor and a memory card port, the QS5 offers repeatable and traceable cleaning performance with consistent cleaning results every time.

Please note: 3 year warranty included when you order with us\*

*\*1 year as standard, 2 years extra warranty applicable if registered online with Ultrawave within 30 days of placing your order*

|                          |  |
|--------------------------|--|
| Overall Dimensions       | 345 x 200 x 265mm (WxDxH)  |
| Internal Tank Dimensions | 300 x 150 x 150mm (WxDxH)  |
| Basket Dimensions        | 275 x 125 x 90mm (WxDxH)   |
| Capacity                 | 4.5Litres  |
| Constrution              | 304 Stainless Steel Tank   |
| Control                  | Digital Control<br>Time (0 - 99minutes)<br>Temperature (Ambient to 70°C) |
| Operating Feaxtures      | Memory Card Port For Digital Traceability<br>Water Level Sensor          |
| Ultrasonic Power         | 100W   |

|                     |                    |
|---------------------|--------------------|
| Operating Frequency | 32 to 38kHz        |
| Heating Power1      | 50W                |
| Required Services   | 13amp 3 Pin Socket |

## ADDITIONAL INFORMATION

|                                 |           |
|---------------------------------|-----------|
| <b>Animal</b>                   | Mixed     |
| <b>SKU</b>                      | 116-617   |
| <b>Brand</b>                    | Ultrawave |
| <b>Ultrasonic Bath Controls</b> | Digital   |
| <b>Ultrasonic Bath Litres</b>   | 12.5      |

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights