

## SET OF 10 MIXED ELECTRODES

GIMA-30500



Price: £106.66

Date Printed: 04/10/23

Category: Accessories > Electrosurgery Units

### PRODUCT DESCRIPTION

## Set of 10 mixed Electrodes for GIMA Electrosurgery Unit

#### Includes:

- Electrode wire 0.2 - straight (code 30501) - Reference number: 1
- Electrode wire 0.2 - angled 45° (code 30502) - Reference number: 2
- Electrode bend - Ø 4 - straight (code 30503) - Reference number: 3
- Electrode bend - Ø 8 - straight (code 30504) - Reference number: 4
- Electrode hook - angled 45° (code 30505) - Reference number: 5
- Electrode wire 1 mm.- angled 45° (code 30506) - Reference number: 6
- Electrode straight slip-knot - straight (code 30507) - Reference number: 7
- Electrode angled slip-knot - straight (code 30508) - Reference number: 8
- Electrode ball point - Ø 3 mm - angled 45° (code 30509) - Reference number: 9
- Electrode ball point - Ø 3 mm - straight (code 30510) - Reference number: 10

### ADDITIONAL INFORMATION

**Animal**

Mixed

**SKU**

GIMA-30500

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights