

**POUR FILLER ADAPTOR**

102-931-1

**PRODUCT DESCRIPTION****Pour Filler Adaptor For Funnel / Pour Fill Style Vaporisers**

Fit a pour fill adaptor to your bottles of volatile agent to prevent spills and inhalation. A simple device to help keep you within the COSHH limits.

The adaptor facilitates the filling of screw /pour fill / well fill vaporisers. Helps transfer agent from the bottle into the vaporiser, minimising spillage/drips etc. Adaptors are agent-specific colour coded for the appropriate anaesthetic agent.

- Purple = Isoflurane
- Red = Halothane
- Yellow = Sevoflurane

**ADDITIONAL INFORMATION**

<b>Animal</b>	Small Animal
<b>SKU</b>	102-931-1
<b>Fill Type</b>	Funnel Fill

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights