

EQUINE BLACK RUBBER REBREATHING BAGS

104-230-BUNDLE



Category: Anaesthesia Consumables > Anaesthesia Consumables > Anaesthesia Accessories > Anaesthetic Circuits > Anaesthesia

Price: £241.10 - £266.50

Date Printed: 17/04/24

PRODUCT DESCRIPTION

Ideal for manual patient ventilation or in connection with anaesthesia apparatus

- Equine Specific
- 50mm Thick neck
- Manufactured from black antistatic natural latex rubber
- Complete with tail and ring (excludes open ended version)
- Replace when surface has because sticky, ripped or brittle

PRODUCT OPTIONS & ADD ONS

Bag Size

1x 10-15 Litre - 50mm Neck

1x 30 Litre - 50mm Neck

ADDITIONAL INFORMATION

Animal	Equine
SKU	104-230-bundle
Catagory Sub Heading	Range of Sizes Available
MRI Compatible	No

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights