

## RESUSCITATION BAGS

104-199\*-BUNDLE



## PRODUCT DESCRIPTION

### For temporary or emergency ventilation of patients.

Designed for ease of use, the intersurgical resuscitator is available in 2 different sizes. The high-quality, robust construction provides the reliability needed to withstand use in all scenarios.

Supplied as a fully assembled, ready-to-use product.

#### Key Features:

- Handle - allowing one-handed ventilation
- Internal diameter of connection will fit ET tube or face mask
- Add on an Oxygen Supply Kit for use within kennels
- Low inspiratory and expiratory resistance
- Highly efficient intake valve minimises CO<sub>2</sub> retention
- Valved diverter cap to minimise dilution of oxygen concentration
- Pressure relief valve - 40cm H<sub>2</sub>O on selected bags
- Slip-free texture and arrow-shape vinyl bag with foldable design
- Clear mask for patient visibility

## ADDITIONAL INFORMATION

Animal

Small Animal

**SKU**  
**Category Sub Heading**  
**MRI Compatible**

104-199\*-bundle

Range of Sizes Available

No

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights