

**BURTONS HD VIDEOSCOPE 13MMOD X 3500WL**

END-KW1335

**PRODUCT DESCRIPTION****Burtons HD Videoscope 13mmOD x 3500mmWL - ideal for mixed practice**

Ideal for small animal upper and lower Gi tract and some URT in medium to large breed dogs, aquariums, zoos, URT in equine, guttural pouches, and cystocopy

**Key Features:**

- 13mm outer diameter
- 3500mm working length
- 2.5 mm biopsy channel
- 4-way angulations 210 degrees up / 100 degrees down
- 90 degree left/right

**Supplied in transportation case with:**

- Biopsy forceps x 2
- Cleaning brush x 2
- Flushing adapters
- Water bottle
- Spare piston kit
- Spare biopsy cap
- Waterproof cap 3 x 2 viewing and recording screen
- Video cables SD card
- User Manual

**ADDITIONAL INFORMATION**

<b>SKU</b>	END-KW1335
<b>Brand</b>	Burtons
<b>Animal</b>	Mixed

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights