

SUNTECH VET20 SPOT CHECK BLOOD PRESSURE MONITOR

SUN-VET20-CONFIG*



BURTONS SunTech Range		SUNTECH VET20 Spot Check BP Monitor	SUNTECH VET25 Spot Check BP Monitor + Pulse Ox + Temperature	SUNTECH VET30 BP Monitor + Pulse Ox + Temperature
MOTION TOLERANT MEASUREMENTS	⊙	⊙	⊙	⊙
SELECTABLE BP AVERAGING	⊙	⊙	⊙	⊙
EDITABLE CLINICAL ALARMS	⊙	⊙	⊙	⊙
SMALL ANIMAL & EQUINE MODELS	⊙	⊙	⊙	⊙
AVAILABLE IN 3 X COLOURS	⊙	⊙	⊙	⊙
BLUETOOTH FOR DATA TRANSFER	⊙	⊙	⊙	⊙
PC APPLICATION INCLUDED	⊙	⊙	⊙	⊙
BATTERY TYPE	AA ALKALINE	LITHIUM ION	LITHIUM ION	LITHIUM ION
BP SPOT CHECKING	⊙	⊙	⊙	⊙
CONTINUOUS BP MONITORING	⊙	⊙	⊙	⊙
TEMPERATURE MEASUREMENT	⊙	⊙	⊙	⊙
SPO2 MEASUREMENT	⊙	⊙	⊙	⊙
Find out more at burtonsveterinary.com		was £1380 now just £995	JUST £695	JUST £795



PRODUCT DESCRIPTION

Suntech blood pressure monitor for spot checking companion animals

The portable SunTech® Vet20 quickly and easily provides reliable blood pressure (BP) measurements on awake companion animals, ensuring a simpler assessment process for veterinarians and their staff members.

Utilizing SunTech's trusted Advantage VET BP technology, the Vet20 quietly performs motion-tolerant, oscillometric BP measurements, delivering accurate results. The intuitive touch screen interface includes all the necessary features for reviewing and analysing BP data.

Alarms can be silenced to avoid frightening the animal and saved measurement data is always retained when the device is powered off.

Watch our Product Spotlight video on the Suntech VET BP Monitor [here](#)

Key Features

- Reliability - As a result of SunTech's motion-tolerant technology, the Vet20 has a high success rate for BP readings.
- Simple to use - Intuitive, touch-screen interface has all necessary features for reviewing and analysing BP data.
- Accurate - BP algorithm was specifically developed for companion animals by experts at a leading vet school.
- Silent - Quietly measures BP and alarms may be silenced to avoid frightening your patients.
- Fast - Take multiple measurements in a matter of minutes!
- Portable - Battery-powered device can be easily moved between exam rooms.
- Memory - The Vet20's advanced memory ensures data is always retained when the device is powered off.
- Averaging: Users can select specific measurements from the device memory and perform automated BP averaging.

Complete package includes:

- Vet20 system monitor
- Protective armour in Flamingo Pink, Peacock Blue, or Tree Frog Green
- Veterinary blood pressure cuffs sizes 1-6
- Patient hose, 4ft (1.2m)
- AC adapter

Specifications:

- Systolic: 40 – 265 mmHg
- Diastolic: 20 – 200 mmHg
- Mean Arterial Pressure: 27 – 222 mmHg.
- Heart Rate: 25 – 300 BPM
- Device Weight: 1.3lbs (with battery), 590 grams (with battery)
- Dimensions: 6.25" x 5" x 5.25", 15.9cm x 12.7cm x 13.3cm
- Warranty: 2 years on monitor, 90 days on cuffs/hose

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	SUN-VET20-config*
Category Sub Heading	Select your colour
Monitor Parameters	IBP, NIBP

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights