

**RED RUBBER E.T. TUBES**

104-01\*CONFIG

**PRODUCT DESCRIPTION****Quality sterilisable Magill endotracheal tubes complete with balloon****Key features:**

- Cuffed and uncuffed option
- Plain design
- Magill style
- Autoclavable to 134 °c
- Reusable
- Do not dry heat

**ADDITIONAL INFORMATION**

<b>Animal</b>	Small Animal
<b>SKU</b>	104-01*config
<b>Catagory Sub Heading</b>	Range of Sizes Available
<b>Tube Material</b>	Rubber

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights