

600 GIVING SETS + FREE BURTONS IVY INFUSION PUMP

129-400-PROMO-600



PRODUCT DESCRIPTION

Easy to use, compact infusion pump, perfect for your practice!

Promotion includes:

- 600 Super Strength Giving Sets
- Ivy Infusion Pump

The IVY vet infusion pump from Burton's is the perfect addition to your practice. Its large, easy to your interface allows for simple set up and monitoring. The compact and lightweight design makes sure the IVY remain unobtrusive when used across the practice.

The Ivy Veterinary infusion pump is part of our Infusion Perfection range, a selection of affordable, veterinary specific, high-quality products and consumables. The range is fully supported and can be calibrated by Burton's in-house service department

Features & Functions

- **Compact design** - The compact and lightweight design saves space but also incorporates a Large & colourful LCD display
- **Added safety** - The built-in Anti-reverse rotation detection and flow stop feature helps reduce the risk of accidental free flow prevent free flow when pump door is open.
- **Anti-bolus function** - Diminishes the volume of unwanted Bolus injected to the patient after removal of the occlusion cause.
- **Panel Lock** - The panel lock feature helps prevent unauthorised changes to any instrument setting

User-friendly operation

Soft key design, easy to operate
Load directly the last infusion rate and volume limit
Adjustable buzzer volume (3 levels)
Top indicator light reflects the operation status as well as alarm case

Night mode (key tone off screen dimmed, top indicator light off) to ensure good rest of removable pump body for easy cleaning

Versatile functions

2000 history records - Real time display

90° rotation pole clamp convenient for horizontal bar, vertical IV pole & platform use

ADDITIONAL INFORMATION

SKU	129-400-PROMO-600
Animal	Small Animal

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights